


Advancements in nanopore sequencing allow in-depth characterization of rAAV vector batches comparable to SMRT™ sequencing

K Breunig, F Dunker-Seidler, F Sonntag, M Hoerer, RC Feiner

Abstract
Long-read sequencing is highly suited for thorough analysis of encapsulated DNA, which is crucial for optimizing gene therapy platforms especially with regards to safety. We utilized nanopore sequencing (Oxford Nanopore Technologies) to profile a rAAV9 batch produced in a 50-liter bioreactor using our proprietary EpyQ™ plasmid system. We evaluated three methods for converting single-stranded to double-stranded DNA, all resulted in consistent DNA impurity distributions. Additionally, we compared recent nanopore advancements, including V14 chemistry and dorado basecalling, with the established V9 chemistry and guppy basecalling, observing significant improvements in





read quality and alignment accuracy. Our nanopore data showed high vector quality (> 95% vector mapping) with minimal host cell DNA (HCD) impurities. Comparison with SMRT™ sequencing (PacBio) confirmed similar HCD and helper plasmid impurity levels, though SMRT™ sequencing detected more plasmid backbone reads, possibly due to differences in chimeric read assignment. Truncation analysis identified consistent vector hotspots, in agreement with an independent outsourced analysis based on the SMRT™ data.



Conversion of ssDNA to dsDNA with

- Annealing
- Reverse Transcriptase
- DNA Polymerase I

For comparison of conversion methods see ASGCT 2023 #1428



Summary

In summary both sequencing methods provide very similar data that prove the high vector quality of the EpyQ™ manufacturing system. Our research highlights the advantages of nanopore sequencing for the comprehensive analysis of DNA encapsulated in AAV particles. Profiting from recent technological advancements, as well as our internal bioinformatics pipeline, we produced data comparable to SMRT sequencing.

However, nanopore sequencing bears substantial benefits, especially during early development as it requires far lower DNA input and has a very quick turn around time. This makes nanopore sequencing an optimal choice for in-depth characterization of AAV-packaged DNA, combining efficiency and precision in a streamlined workflow.

