

Streamlined pre-formulation screening with minimal sample requirements & 2-day turnaround to enable high efficiency AAV formulation development

Philipp Beck, Astrid Stockinger, Sonya M. Schermann

Abstract


Material constraints are a limiting factor in early stage AAV gene therapy development programs. This is especially true for initial, broad formulation studies, that consume substantial amounts of material and time. We have established a low-volume, high-throughput pre-formulation screening approach using a simple proprietary buffer matrix that spans a broad pH range (5 – 8) and ionic strength levels (using 100 – 310 mM NaCl). We show that rebuffing of AAV (exemplary serotype AAV3B/8 hybrid) using a commercially available 96-well spin-plate loaded with size-exclusion resin provides a rapid and reliable means of rebuffing. With only 200 µL per condition, the material from this rebuffing campaign was used to generate thermal (nanoDSF) and colloidal (DLS) stability data along with assessing subvisible/visible particle content (BMI), DNA ejection (using SYBR GOLD) and titer recovery (Stunner).

Since all methods are compatible with the 96-well format, results from this screening approach are available within hours of the rebuffing procedure with minimal manual handling. For the exemplary serotype used in this study, a buffer at pH 8 (20 mM TRIS, >190 mM NaCl) initially provided satisfactory results, which is in line with published data. However, DNA ejection progressed steadily at the tested temperature (37°C) over 2 days compared to all other tested buffers, which exhibited a much less pronounced payload release. Thus, development of a liquid formulation for the serotype in question should focus on buffers at lower pH levels to mitigate loss of encapsulated DNA during storage. Overall, the wealth of data generated from this method with just 6.5 mL of sample AAV (>2E12 vg/mL) can provide meaningful insights towards buffer and ionic strength preferences of a given serotypes and further inform rational formulation development. The approach is capsid agnostic and can be applied to any AAV serotype or other modalities based on protein nanoparticles.

Colloidal stability

Dynamic Light Scattering


Low-salt conditions lead to aggregation after rebuffing, which is evident as an increase in dispersity of the sample (yellow fields). pH levels of 7 or 8 resulted in less polydispersity after rebuffing compared to pH 5 or 6, potentially hinting at better resistance against aggregation during downstream processing and storage.



Subvisible particles

Backgrounded Membrane Imaging

Evaluation of subvisible particle content is usually constrained to high-volume formulation studies at later stages due to high volume requirements of e.g. light obscuration. Backgrounded membrane imaging requires only microliters of sample and provides an insight into size & shape of particles. High-salt conditions lead to reduced particle formation due to rebuffing across pH levels (blue fields).



Identical 384-well plate used for evaluation of AAV aggregation upon heat-stress & simultaneous DLS assessment

Heat stress resistance: 20°C → 90°C over 2 h

Monitoring temperature-induced onset of AAV aggregation (T_{onset}) showed that ionic strength in addition to pH are main contributors to colloidal stability upon heat stress (right). This result confirms published data[1], validates our platform, and is in line with nanoDSF data (shown below).

pH	100	130	160	190	220	250	280	310
8	31.4	38.4	27.2	59.0	64.5	62.8	70.6	71.1
7	30.5	49.4	64.3	61.4	61.3	63.5	62.8	62.6
6	46.6	46.7	38.5	53.7	69.2	64.5	83.7	80.3
5	49.5	46.5	47.7	49.9	54.6	61.5	62.9	64.2

Your AAV
+ Rebuffing spin plate = Results within 2 days

Only 7 mL of > 2E12 vg/mL required
Exemplary buffer set: 20 mM TRIS/Phosphate/Histidine/Acetate


- 1) Equilibrate spin-plate with buffers
- 2) Apply AAV, spin down using centrifuge
- 3) Distribute rebuffed AAVs using multi-tip pipette

60 µL, 40 µL, 30 µL, 10 µL

DNA-Ejection

SYBR-Gold Fluorescence


Measurement of SYBR-Gold fluorescence as an indicator of ejected DNA at 37 °C over 24 hours is performed in simple thermocyclers – no need for specialized equipment. For all buffer levels, rapid DNA ejection is observed in the first 10 hours of incubation. Subsequently, there is only minor additional payload release evident at that temperature except for the condition at pH 8, which exhibits a steady increase in SYBR-Gold fluorescence. Overall, the extent of DNA ejection is dependent on NaCl concentration in each buffer level tested, with higher salt concentrations leading to less ejection.



Capsid stability

nanoDSF


Capsid unfolding temperatures are pH-dependent and not influenced by ionic strength. A pH level of 6 resulted in the highest capsid stability during heat stress. Formulations at pH 7 exhibited the lowest resistance to heat with melting temperatures at 76°C, irrespective of ionic strength. In line with previously reported data, a drop in pH from 7 to 6 leads to a marked increase in capsid stability of 10°C.



Recovery after rebuffing

DLS/UV-Vis using Stunner

Assessment of capsid recovery after rebuffing using Stunner's DLS/UV-Vis combination revealed a strong dependency on ionic strength at pH 8. Downstream processing should thus avoid "salt valleys" especially at this pH level – our DSP-poster P0031 shows how we applied this knowledge to get to >50% vector-recovery in the purification process. NaCl-levels >190 mM did not further enhance the recovery at pH 5/6/7.



Summary

The poster outlines how a wealth of formulation data can be generated within only 2 days from minimal amounts of sample. Using a rebuffing approach via 96 well spin-plates enables sample handling with either a pipetting robot or (as in our case) or a multi-channel pipette. Compared to other methods of rebuffing (Slidealyzer

cassettes/cups, spin-filters "amicons", PD-10 columns), the presented workflow is fast and tailored to the low-volume 96-well-compatible methods that are already in place in most laboratories. Overall, the workflow is capsid agnostic and can be applied to any AAV serotype or other modalities based on protein nanoparticles.



Ascend Advanced Therapies is a global gene to GMP development partner

Germany | United Kingdom | United States

www.ascend-adv.com
Philipp.beck@ascend-adv.com

Aim higher